

BI Sobrietor®

If the Problem is Alcohol,
the Sanction is Sobrietor



- New variable security levels reduce false alerts and minimize attempts to fool the system.
- To make remote testing consistently reliable, voice prompts are incorporated into Sobrietor that provide the offender with clear, unambiguous instructions.

More than 1.5 million DUI arrests are made¹ and approximately three million alcohol-related violent crimes are committed each year.² Now it is possible to sentence the individuals who commit these crimes to sobriety with the newest multi-language alcohol monitoring technology from BI Incorporated. With Sobrietor®, an offender can be required to remain sober. Since monitoring is performed remotely using a telephone connection, compliance with the alcohol restriction is verified each and every day.

Easy Enrollment

The Officer Control Unit makes offender enrollment simple and fast. The unit's visual display and keypad guide the officer step-by-step through client enrollment, minimizing training. Plus, the Officer Control Unit allows the officer to use Sobrietor to verify voice enrollment and perform on-demand field sobriety tests.

Sobrietor is available for English and Spanish speaking offenders. Multi-language units are available that have this capability, and the officer determines the language setting through the Officer Control Unit.

Simple, Reliable Testing

To make remote testing consistently reliable, voice prompts are incorporated into Sobrietor that provide the offender with clear unambiguous instructions. For example, if the individual is speaking too softly, Sobrietor prompts the offender to "speak louder." Voice prompts ensure that the offender is instructed through each step of the test.

Officer-controlled Security

New variable security levels give the officer the ability to adjust the voice verification threshold for each offender. This feature reduces false alerts and minimizes attempts to fool the system. In addition, Sobrietor stores the last failed voice recognition test, helping to identify uncooperative offenders.

Key Features

BI Sobriector®

Environment

- Temperature: 59-104°F
- Humidity: 15-95% non-condensing

Compatibility

- GuardCenter®
- GuardWare®
- Host Products
- BI 9000™ Field Equipment
- BI HomeGuard® 200
- BI VoiceID®

Battery Back-up

- Up to 12 hours
- Up to 4 tests

Fuel Cell

- Calibration recommended every 6 months

Language


- English
- Spanish



innovation

expertise

BI builds and operates an integrated suite of reentry services and technology-based monitoring solutions. We are focused on helping corrections agencies break the cycle of recidivism and improve community public safety.



Scope: The design, production, installation and support of EM services and products in the Boulder, CO and Anderson, IN facilities.

Accurate Testing

Sobriector utilizes fuel cell technology, the same alcohol consumption monitoring system used in the most advanced Portable Breath Test devices. Because Sobriector measures the Breath Alcohol Content (BrAC) from the breath test, the offender's intoxication level is reported to the officer.

Sobriector and Radio Frequency Monitoring: A Powerful Combination

Sobriector may be used as a stand-alone sanction or in conjunction with Radio Frequency Monitoring. In stand-alone mode, Sobriector tests the offender randomly throughout the test period, as specified by the officer. When used with Radio Frequency Monitoring equipment from BI, the offender may be tested immediately upon entering the house, and randomly thereafter. In either case, only one phone line and one power outlet are required.

Monitoring of Sobriector tests is performed at the BI GuardCenter®. Trained BI personnel perform testing procedures specified by the agency, and report all results. Agencies may also utilize an on-site host system allowing them to control compliance monitoring internally. BI works closely with every customer to determine which monitoring system best meets their needs.

Sobriector has the ability to help individuals break the relationship between alcohol and criminal behavior. Many crimes can actually be prevented by requiring mandatory sobriety. BI offers remote alcohol testing as an essential component of our continuum of graduated sanctions. Sobriector is an effective method to monitor an offender's compliance and progress, and can be an important part of BI's integrated supervision and treatment products and services.

For more information about BI Sobriector call 800-701-5171 or visit www.bi.com.

1. Bureau of Justice Statistics, Alcohol and Crime: Analysis of National Data on the Prevalence of Alcohol Involvement in Crime, 1998.

2. Federal Bureau of Investigation, Uniform Crime Report, 2000.

Sobriector®, GuardCenter®, BI HomeGuard® 200, BI 9000™, GuardWare®, and BI VoiceID® are Trademarks of BI Incorporated. 07/05